
MAXXLINE EKSPECTORATSET

EXPECTORATED SPUTUM CONTAINER

STERILE INNER CONTAINER

SWAB



Product Data

General laboratory product for collection and evaluation of sputum samples.





	Inner container		Outer container		Swab
Ref. no.	07.061.300		08.212.1003		1131899
Description	Inner container for sputum Cup with screw cap		Outer container for sputum With snap cap		Cellulose swabs
Dimension (mm)	Ø65 / H 22,8		Ø72,8 / H 26,7		40 x 50
Volume	24 ml		N/A		12 ply 21-17 g/m ² per ply
Material	Cup	Cap	Cup	Cap	Cellulose Squeezed, crimped, pinched at the edges. Straight-edged twisted on a spindle
	Polystyrene (PS)	Polyethylene High-Density (HDPE)	Polyethylene High-Impact (HIPS)	Polyethylene Low-Density (LDPE)	
	Physiologically safe, biologically inert, free from heavy metal				-
Color	Cup	Cap	Cup	Cap	-
	Clear	Natural	Natural	Yellow	Oxygen bleached
Sterilization	Sterile		Non sterile		Non sterile
Temperature	-20°C to 60°C		-10°C to 60°C		-
Quality Control	Testing of attributive and variable characteristics in accordance with the valid specification				-
Storage	To be stored in closed original packaging Dry, at room temperature				
Expiration time	5 years				
Autoclavable	No				
Regulations	Directive 98/79/EC		DIN EN ISO 13485		DIN EN 13485
Other information	Approved by Danish Technological Institute after ADR2011 and P650. For single use only.				

Shipping info

Catalog no.	AL-450002-S
Box packaging	50 bags x 5 pcs
Box weight (kg)	10,000
Box dimensions (WxDxH) (mm)	400 x 600 x 385
Eu pallet packaging	20 boxes (1 lay = 4 boxes)
EU pallet weight (kg)	200,00
EU pallet dimensions (WxDxH) (mm)	800 x 1200 x 2040

User manual

User manual version 4.2.2 is placed on the top of the inner container.

Brugermanual Usermanual	<p>1. Tag det inderste bæger ud, og skru det hvide låg af. Take out the inner container and screw off the white cap.</p> 	<p>2. Placér ophostet prøve. Place the sample.</p> 	<p>3. Skrú låget stramt på igen. Screw the cap on tightly.</p> 	<p>4. Kom det lille bæger tilbage i det store bæger. Place the small container in the big container.</p> 
--	---	---	--	---

+45 76307600 • info@dacos.dk • www.dacos.dk

